

torq.control®



**Universal
torque wrench.**

ANTHOGYR

A Straumann Group Brand



torq.control®

For controlled screwing.

Torq Control® is the Anthogyr universal torque wrench ensuring precise tightening, crucial for implant success, regardless of connection type or access constraints. Designed for all prosthetic procedures, it offers optimal safety with automatic declutching.

A must-have, especially for full arch restorations.



ANTHOGYR 
A Straumann Group Brand

torq.control®

HIGH PERFORMANCE.

Experience versatility and precision with a solution designed for every dental implant system. Combining proven technology with ease of use, Torq Control® ensures optimal torque management, effortless maintenance, and maximum safety.

COMPATIBLE WITH ALL IMPLANT SYSTEMS!


- Universal connection for all dental implant systems
- 7 torque values, from 10 to 35 N.cm, adapted for all protocols

DESIGNED FOR EFFICIENCY

- Micro-head, 100° angulated for easy access and optimal visibility, even in posterior sectors



Universal torque wrench.



ERGONOMICS AND EASE OF USE.

MAXIMUM COMFORT OF USE

- Manipulation 100% outside of mouth
- No lateral movement

SO EASY TO HANDLE


Only 130 g for a better freedom of movement

CONVENIENT MAINTENANCE

- One-piece design
- Sterilizable in autoclave, 135°
- Repairable

DID YOU KNOW?

Trusted by professionals since 1997, Torq Control® remains a must for its precision and ergonomics. Still one of Anthogyr hero products!



Watch Torq Control® in action

RELIABILITY AND SAFETY AT THE CORE.

PERFECT PRECISION IN ALL SAFETY

- High precision torque tightening
- Automatic declutching when the selected torque is reached

GUARANTEE OF QUALITY

- Made in the Mont-Blanc Valley, France, since 1997
- High quality manufacturing



MADE IN FRANCE
An iconic product since 1997!



campus

byANTHOGYR

IMMERSION 360

WELCOME TO ANTHOGYR.

STAY TUNED!



Step into our world and explore a 360° view of our expertise and operations at our Sallanches site in the stunning Chamonix Mont-Blanc Valley.

A unique opportunity to meet our teams, uncover the behind-the-scenes of our industrial site, and dive into the secrets of Anthogyr's cutting-edge implantology solutions—all set against a breathtaking backdrop.

FOLLOW US.



www.anthogyr.com

ANTHOGYR

A Straumann Group Brand

